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---

### EN Seat cushion for pelvic floor training – Instructions for use

**Dear Customer,**

Thank you for choosing our seat cushion!

The Instructions for Use should help you to get the best use out of your product. Please read the instruction manual carefully – particularly the usage and safety instructions – before using this product for the first time.

We will be glad to answer any questions and to receive any suggestions you might have!

*Your ELANEE Team*

Please observe and keep the instructions for use and the packaging!

## **1. Field of application**

These instructions for use are applicable for the following product:

**REF** 708-V1 seat cushion for pelvic floor training

## **2. Product content**

- 1 seat cushion for pelvic floor training
- 1 hand pump
- 1 Instructions for use including an exercise programm

## **3. Purpose and intended use**

The seat cushion increases comfort when sitting and at the same time provides a supportive function for training the pelvic floor.

Thanks to its ergonomic shape, the cushion enables an active and conscious sitting posture, which helps to strengthen and activate the pelvic floor muscles.

The seat cushion also promotes awareness of the pelvic floor muscles, which can lead to improved posture and an increased sense of well-being in everyday life.

## **4. Areas of application:**

The seat cushion is suitable for passive and active training in the pelvic floor, stomach and back muscle regions. This helps to stabilise the spine and pelvic floor, while also training coordination and balance.

## **5. Functionality**

The seat cushion's operating principle is that of "active sitting". The cushion is filled with air and when you sit on it, it creates minor instability in the pelvis. The pelvic floor muscles counteract this instability and are thereby permanently in use and stimulated while sitting. This constant movement trains the pelvic floor muscles.

The seat cushion for pelvic floor training is the perfect addition to pelvic floor exercises or training with pelvic floor training aids (e.g. ELANEE phase I or phase II active pelvic floor training aids).

## **6. Usage**

1. Check that the cushion is in good condition.
2. If there is not enough air in the cushion, pump it up using the enclosed hand pump. To do so, screw on the pin of the pump and insert this into the valve on the cushion. Please take care not to damage the sitting cushion with the

pin! Pump air into the cushion until it is full enough. To remove air from the cushion, simply insert the pin into the valve without the pump and press gently on the cushion until enough air has been removed.

3. Place the cushion on a firm, flat sitting surface with the knobbly side facing downwards. Unpadded surfaces are particularly suitable.
4. Sit on the cushion.
5. Stay sitting on the cushion for as long as it feels comfortable and/or conduct 10 minutes of exercises every day.

**Tip:** You can find exercise ideas for effective pelvic floor training using the sitting cushion in the exercise programme (in DE, EN, FR) which can be found in this booklet.

## 7. Safety information

Only use the seat cushion on solid and level seating!

Avoid contact with sharp or pointed objects in order to prevent damaging the cushion!

Do not use abrasive cleaning agents or disinfectants containing alcohol or chlorine and do not use hot air for drying!

## 8. Safekeeping, service life, storage and transport

- The sitting cushion has been conceived exclusively for use in dry indoor areas. It should also be stored in a dry place.
- The sitting cushion is intended for repeated and long-term use.
- We recommend replacing the sitting cushion after a maximum of five years, if used normally.

## 9. Cleaning and disposal instructions

- The seat cushion can be cleaned with a mild cleaning agent and a moist cloth.
- Please separate the individual components of the packaging according to different materials and recycle them appropriately. Please observe the local and regional disposal regulations in the process.
- The sitting cushion can be disposed of with normal household waste.

## 10. Material composition

100% polyvinyl chloride

## 11. Explanations of the symbols used

Manufacturer

Observe the instructions for use!

Batch number

Article number

Attention!

Labelling of the packaging material

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## Exercises for professionals


### Exercise 5 - Foot spin

**Starting position:** Place the cushion on the floor. Place yourself on the cushion with both your feet. This exercise is similar to exercise 3 - however requires more balancing sense.

- ⚠ Be careful when sitting yourself on the cushion. Find a safe holding-on point as support such as the back of a chair or a door frame!
- Let the feet rotate in the clockwise direction - forward, to the right side, backward and to the left side. There should be a circular movement. Try to keep your body as stable as possible and shift only your centre of gravity.
  - Do ten circles clockwise, then at least 1 minute of rest followed by 10 circles counter clockwise. Repeat the entire sequence three times.

### Exercise 6 - Bridge

*Starting position:*

**Starting position:** Place the cushion on the floor. Sit on the cushion. Support yourself behind your spine with both the arms. Your feet will remain together on the floor.

- Lift the legs approx. 20 cm above the floor. Keep the spine straight and try to keep the body stable.
- Hold this position for approx. five seconds.
- Move back to the starting position. Hold the starting position for five seconds.
- Repeat the exercise ten times.

## Exercises for experts

### Exercise 3 - Hip rotation

**Starting position:** Place the cushion on a stool. Sit on the cushion. This exercise is a combination of exercise 1 and exercise 2.

- Let the pelvis rotate in the clockwise direction - forward, to the right side, backward and to the left side. There should be a circular movement.
- Do ten circles clockwise, then at least 1 minute of rest followed by 10 circles counter clockwise. Repeat the entire sequence three times.

### Exercise 4 - Balance


**Starting position:** Place the cushion on the floor. Sit on the cushion. Bend your legs, your feet will remain together.

- Lift your legs. Keep the spine straight and try to keep the body stable.
- Bring your hands towards the knees. Hold this position for approx. five seconds.
- Move back to the starting position. Hold the starting position for five seconds. Repeat the entire sequence ten times.

## EN Seat cushion for pelvic floor training – Exercise program

### Exercises for beginners


#### Exercise 1 - Cat's arch

**Starting position:** Place the cushion on a stool. Sit on the cushion.

- Push your pelvis forward. Your back will have a slight hollow.
- Push your pelvis backward. Your back will be in a cat's arch.
- Repeat the exercise ten times. Take a break for at least 1 minute and repeat the exercise. Repeat the entire sequence three times.

#### Exercise 2 - Hip swing

**Starting position:** Place the cushion on a stool. Sit on the cushion.

- Push the pelvis to the left. The right hip will go up.
- Push the pelvis to the right. The left hip will go up.
- Repeat the exercise ten times. Take a break for at least 1 minute and repeat the exercise. Repeat the entire sequence three times.

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