Beckenboden- Trainingshilfen Phase I aktiv

- ⒹE Übungsheft S. 02
- ⒺN Leaflet of Exercises ..p. 20
- ⒻR Livret d'exercicesp. 36

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EN Pelvic Floor Training Aids Phase I active – Leaflet of Exercises

Preface

Welcome to the Training Course!

You have decided to improve your quality of life. A strong pelvic floor not only counteracts a weak bladder and a dropping of the bladder but also supports the female basis. The strength coming from your core will straighten your internal alignment, protect your back and will positively influence your body posture altogether as well as your attitude towards life. By practicing this holistic pelvic floor training you lay the foundation for active long-term health care. And you will also gain more energy and strength daily.

Even by only wearing the ELANEE Pelvic Floor Training Aids Phase I active - as described in the instructions - you will be much more aware of your pelvic floor and strengthen your muscles. If you combine this with the training program presented here, you will benefit even more. You may also discover that it is also a very pleasant exercise.

The pelvic floor consists of three layers of muscles. By using the cones you will primarily exercise the first layer - sphincter - and the second layer that encloses the vagina. This is particularly stressed by giving birth and must be strengthened sustainably.

By doing these exercises you will additionally train the third pelvic floor layer effectively. This layer protects the uterus, intestines and bladder, and plays an important part in keeping an upright posture. Once you build up this part of the pelvic floor then you will be stable from the body basis up. You will not exhaust yourself so easily and your mo-

vements will become more harmonious and more vibrant at the same time.

Structure of the Training Program

As a pelvic floor therapist I meet women every day who find it difficult to become aware of the „secret“ muscular system. The vaginal cones may be of great help in this respect. Apart from the training effect that you will achieve with these cones, they are also wonderfully suitable to make you aware of your pelvic floor muscles. You will feel exactly where you have to contract muscles and you will be able to train more precisely every time.

The training program starts with two awareness exercises, followed by two highly effective exercises to build up strength and finally two exercises that will help you to integrate an active pelvic floor into your normal every day movements.

Training Recommendation

If you start with exercises of 10-15 minutes daily and wear the cones for a further 10-15 minutes after a few hours in between, you should be able to notice a distinct improvement after a few weeks. By incorporating the exercises into your daily routine you should even be able to reduce the active time of training.

This program has been designed such that you can carry out the exercises of a complete unit simply one after the other. Change into comfortable clothes and have everything ready that you need. Use the same cone for the exercises that you wear at present according to the instructions for use.

And here we go:

1. Awareness

The better you can feel your pelvic floor the better you will carry out the exercise - because you activate the correct muscles. Often women contract their abdominal, bottom or thigh muscles intensely, which does not help the pelvic floor at all. Contracting the pelvic floor correctly and intensely, is still a comparably small movement - and it is enough.


The following exercises will help you to discover your body's most "secret" muscles in detail and get acquainted with them.

Exercise 1: Blinking

Position:

Sit upright and relaxed like a queen (don't slip into a hollow back) at the front half of your chair, spread legs apart at hip-width. You should be able to feel both your sit bones against the chair.

Now position yourself on the cone such that it is placed at the positions outside of your clothing as marked in the graphics. Slide it first beneath the vagina (1). Feel the pressure and blink at the cone several times with a small movement of the muscles. This should be completely invisible at the outside. Does it feel almost like gripping? That is exactly right.

Slide the cone now further towards the back such that it rests beneath the perineal area (2). Try and move the sit bones slightly towards the cone, slightly drawing them together. Again it is correct if the movement cannot be noticed from the outside.

Finally push the cone further still towards the back, until it rests underneath the anus (3) and try again to blink and to grip gently.

Check your buttocks' muscles with your hands. If you can keep them still and despite that something around the cone moves distinctly towards the inside then you are actually activating your pelvic floor quite expertly already.


Exercise 2: Gripping

Position:

Lying comfortably on your back, upper body slightly raised, either with a firmer cushion or on your sofa, provided it is not too soft and has a suitable arm rest. Place your bent legs apart at hip width.

If you wear the cone generally, you should push it towards the inside as far as possible such that you do not feel it anymore, just as shown in the graphics. We have deliberately changed this, however, just for this awareness exercise. Only insert the cone so far that it just disappears and keep hold of the cord. Once you tighten up your vaginal muscles you should be able to feel the cone as it is now positioned exactly in line with this layer of muscles.

If you tighten the muscles further the cone will move with the muscle movement; it may even drop out or is drawn in further; try to enclose and move the cone precisely. Play with it by holding it either tighter or only very lightly. Bottom, thighs and abdomen remain still. Above all take care that you do not use your abdomen muscles; the correct movement is very small and „subtle“.

As soon as you have welcomed your “secret” muscles appropriately, push the cone into the vagina until you cannot feel it anymore and continue with your strength-building exercises.

Ideally you should set aside enough time for the awareness exercises at the beginning. You should be undisturbed, alert and relaxed to give your pelvic floor your full attention. Once you have obtained a good feeling for your strength you can reduce the time for this first training section.

2. Building up Strength

Wearing the pelvic floor training aids encourages the holding and thus a static build-up of muscles. If you carry out the exercises in this leaflet regularly you will gain a pulsing quality by tightening up and gently relaxing the muscles. In doing this you achieve a holistic effect and your muscles are built-up more comprehensively.

During the awareness exercises you have blinked and played, but now you are going to work seriously with the following strength-building exercises. Your breathing will support you. Ideally tighten up the muscles when breathing out and relax them when breathing in. If you do not succeed in using this rhythm, take care that you do not hold your breath and seize up. Furthermore, it is recommended not to tighten up the muscles jerkily but increase the strength steadily and also relax steadily - just like a wave comes and goes.

ATTENTION: PLEASE BE AWARE:

If you feel that the exercise is uncomfortable for you and exhausts you and if you get the feeling that when tightening up everything “presses towards the bottom” you are obviously pressing too much with your abdomen muscles. And this is not good for you.

Try to relax a little, only use half as much force as you are able and try and get strength „from the pelvis upwards“. You are doing the exercise correctly, if you get a feeling that “something draws upwards” that it is strong and at the same time somehow light, and that you feel taller.

Exercise 3: Rolling the Pelvis

Position:

In a slightly raised back position or flat on the floor.

a) Gentle Rolling of the Pelvis:

Keep hold of your trousers in the crotch area. Breathe in deeply and relaxed. When breathing out roll your pelvis upwards: Stretch your back as if it wants to push down towards the heels and pull the pelvis gently towards you by holding on to your trouser material. While breathing in your pelvis rolls back into its initial position. Repeat this rolling movement 10 times in line with your breathing, without tightening up any muscles.

Your back will love this exercise. Always make sure that your back really only stretches out lengthwise! While lying down you will hardly feel the cones.

b) Powerful Rolling of the Pelvis:

Breathe in deeply and relaxed. While breathing out tighten the pelvic floor, roll your pelvis upwards and at the same time press your heels firmly into the floor support. Don't forget to draw the back out lengthwise!

Release the tension while breathing in, and let the pelvis sink back into its initial position.

Repeat 10 times in line with your breathing.

Position:

In a slightly raised back position or flat on the floor as in exercise 3.

Exercise 4: Pressing down the Ball of the Foot

Breathe in deeply and relaxed. Tighten up the pelvic floor when breathing out and press down the front section of your feet. In doing this heels and pelvis will slightly lift up.

It is not necessary that the heels move up (as shown in the picture) - only do what you feel is comfortable for you.

It is, however, important that you do not fall into a hollow back. The buttocks can help a little, avoid moving the abdomen as far as possible. You may perhaps get a feeling that you are firmly embracing the cone with your vagina and sit bones. And this is exactly how it should be.

Release the overall tension while breathing in, and the pelvis sinks back into its initial position.

Repeat 10 times in line with your breathing rhythm.

Tip: Always start and finish your exercises with a soft rolling of the pelvis. Carry out the powerful rolling of the pelvis 10 times and the pressing down of the balls of your feet 10 times. In between, gently roll the pelvis 10 times. The change between powerful tension and gentle rolling helps to build up the muscles.

3. Incorporating the Exercises into your Daily Routine

While you can consider the toning exercises as a type of „dumbbell training“ for the pelvic floor, the following exercises concentrate more on body coordination. You do not have to use a lot of tension, but a good rhythm is more important. The pelvic floor training aids reinforce the training effect.

If you do these exercises regularly you will learn to let your everyday movements such as climbing the stairs for example come increasingly more out of the pelvic floor. In doing this you can exercise even without cones and gain additional elasticity and agility - as the pelvic floor is your body's energy center. If you keep it active by the right kind of movement, you are efficient and cheerful.

Exercise 5: Dance on a Stool

Position:

Sitting down on a stool or a chair with a level seat, not too soft.

Do this 5 times on the right, 5 times on the left. Repeat the whole sequence 3 times. Then 20 times the right and left alternately.

Place a cushion or, better still, a cherry pit cushion on your head. That will help you to remain upright during this exercise and not to move the upper part of your body.

Tighten your pelvic floor while breathing out and slowly lift your right leg just a little off the floor. It would be easier to lift your leg by changing the weight and tilt your upper part of the body to the left. Resist this temptation and use your strength coming from the pelvis to lift the leg. It is helpful to slightly press the left foot onto the floor and turn your left shoulder slightly towards the front.

Try and imagine that you are growing a little every time you tighten the pelvic floor and lift your leg.

Whilst breathing in slowly place the foot back on the ground in a concentrated manner.

Exercise 6: Climbing Stairs

Position:

Stand in front of stairs or a low, safely positioned elevation.

Again place the cherry pit cushion on your head so that you keep yourself upright and do not swing your upper body to and fro. Place your right foot on the stairs.

Tighten the pelvic floor while breathing out, and slowly shift your weight completely onto your right leg until the toes of your left foot only just about touch the floor. Get the strength from the pelvis, even though the thigh is of course needed in this exercise, and again grow a little. While breathing in shift your weight back to the left foot and place it back on the ground.

Do this 5 times on the right, 5 times on the left. Repeat the whole sequence 3 times.


If you do your exercise on a staircase, move up the whole staircase slowly and in a concentrated manner.

4. Relaxing

This is more important than you might think. So that the muscles can be built up after the exercises have finished, they must be supplied with blood and thus also with nutrients.

Indulge yourself!

Relaxation Exercise 1: Enjoyable Stretching

Just do what dogs and cats love to do:

Extensively stretch and luxuriously draw out your arms, legs and back. Let your body decide what's good for it.

Relaxation Exercise 2: Inner Tracing

Cuddle up into your most comfortable position and place your hands on the abdomen. Picture your pelvic floor and inner organs as best as you can. By the way, body awareness reinforces the effectiveness of building-up training of muscles.

Just imagine that you are pleasantly sinking deep into your seat. Hum a deep tone quietly or sigh when breathing out. Allow yourself to get acquainted lovingly with this part of your body. Does it tingle or “buzz” down there? Perfect!

Relaxation Exercise 3: Warming Deep Relaxation

Place a hot water bottle, or a hot cherry pit cushion on your abdomen. There is nothing like letting go deep down - and you have deserved this feeling of wellbeing!

Final Recommendations

Advice on wearing the ELANEE Pelvic Floor Training Aids Phase I active: The more you move with it, the more agreeable it will become for you. While standing or sitting, the pelvic floor muscle system tends to slacken. It tightens up as a reflex when the cones are worn in order to prevent them from gliding out, however, this continuous tension can easily lead to cramp. If women do not feel comfortable with the cones, then this is more often than not the reason. But if you move, wearing the cones will become an inner massage.

Here are some recommendations for your everyday life. Your pelvic floor will always benefit from movement and good posture, whether you are wearing the cones or not.

- If you have to stand for a longer period of time, whether you are ironing or any other activity, it is recommended that you resume a stepping position.
- Avoid a round back when bending down. Simply go in to a stepping position by bending your knee, maintaining a straight upper body.
- Should you have to wait in a queue for a longer period of time, you can shift your weight slowly and invisibly from one leg to the other and your pelvic floor will benefit at the same time.

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GRÜNSPECHT Naturprodukte GmbH