ELANEE

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
Literatur: 1. Egan ME, and Lipski MS: Diagnosis of Vaginitis. American Family Physician 1–14 (September 2000). 2. Sobel JK et al.: Vulvovaginal candidiasis: Epidemiologic, diagnostic and therapeutic considerations. American Journal of Obstetrics and Gynecology V178 (2) (1998).

Instructions for Use - ELANEE Test for vaginal fungal infections

Dear Customer,
Thank you for purchasing our ELANEE Test for vaginal fungal infections. These instructions for use are to help you to use the product properly. Please read the instructions through carefully- in particular the "safety advice" - before using this product for the first time.
We will be glad to answer any questions and receive any suggestions you might have.
Your ELANEE Team
 Please read through and keep the instructions for use and packaging!
Made in Germany

- 1. Scope of Application**
These instructions for use apply to the following product:
REF 725-00 ELANEE Test for vaginal fungal infections
- 2. Contents**
1 test cassette + 1 bag with desiccant,
1 smear swab
- Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia, Italy
- CE 0123** **STERILE R**
- 1 set of instructions for use

3. Purpose and Intended Use
The ELANEE Test for vaginal fungal infections is a quick, simple test which you can do at home yourself and which indicates a yeast (candida) infection in the vaginal secretion.

The ELANEE Test for vaginal fungal infections is suitable for identifying a yeast infection with symptoms such as acute itching, increased discharge which often smells unpleasant or a rash on the outer regions of the vagina.

The presence of a line on the test line indicates a positive result. The presence of an additional line in the control region indicates that the test has been used properly and is working correctly.

4. Area of Application
If you notice one or more of the following symptoms we recommend using the ELANEE Test for vaginal fungal infections:


- Acute itching in the genital area

- Burning sensation in the vagina and outer genital area
- Irritated skin and soreness
- Increased vaginal discharge, often with an unpleasant smell
- Burning sensation in the genital area which may be more severe when urinating

The test is not suitable for any other purposes than the ones named above.

5. How it works
The test for vaginal yeast infections is a point of care test to do at home for the qualitative detection of yeast infections on vaginal smear swabs within 10 - 20 minutes. It is an important advance on improving the diagnosis of women with Vulvovaginal candidiasis (VVC).


- 6. Application**
1. Ensure hygienic conditions (e.g. wash your hands) and have a timer nearby to measure the time.
 2. Remove the device from the sealed pouch and place it on a clean, level surface.
 3. Hold the device with one hand. Carefully remove the aluminium cover of the cap with the other hand.
Note! The cap contains a liquid!
 4. Remove the swab from the wrapper.
 5. Insert the swab one inch (2cm) into the vagina, and rotate for approx. 20 seconds. Pull the swab out carefully!
 6. Dip the swab in the liquid in the cap and rotate it for approx. 20 seconds.
 7. Pull the swab out carefully by squeezing it against the inner wall of the cap. Discard the swab.
 8. Hold the test cassette with one hand. Turn the cap counter-clockwise until it stops (b), and then turn it back to starting point (a).
 9. Repeat step 8 two more times. In the final position the notches on the cap and the white base must be aligned (b).
 10. Read the results after 10 minutes. A line should appear at the control (C) mark. For a positive result, an additional line should appear at the test (T) mark. For a negative or unclear result read again after 10 more minutes..

Positive
Appearance of both lines - detection of a yeast infection.
Note: The intensity of the colour in the test line region (T) should be considered positive even if the coloured line is faint.


Negative
Appearance of only the control line (C). No detection of a yeast infection



Invalid
If both lines do not appear, or if only the test line appear, the test cannot provide a result. It is invalid. A new test should be carried out and the instructions for use must be followed exactly.

7. Safety Advice
 Ensure hygienic conditions (e.g. wash your hands)!

 Please check the packaging before the first use. Only use sealed and undamaged components. If the packaging has been opened or damaged, this may lead to incorrect results.


 Do not use if the sterile packaging is damaged or the use by date has been exceeded! Only use if the packaging is sealed.


 Do not use this test within 7 days of taking medication against vaginal infections.

 Do not use this test within 24 hours of using creams, gels, foams, rinsing solutions or any other products for the vagina.

 Do not use the test during your menstrual period.


 Open the foil packaging directly before use.

 The product is designed to be used once!

 Do not eat, drink or smoke while performing the test.

 Only for use as in vitro diagnostic. Do not consume.

 A negative test result does not mean that there are no other infections. Please see your doctor if your symptoms persist or worsen!


 Side effects may include skin reactions such as reddening of the skin, swelling or burning. If you know that you are sensitive to any of the product components mentioned, we advise against using this product.

8. Storage and Shelf Life
Please store the test for vaginal fungal infections in a dry place at room temperature of 2°C to 30°C, keeping the packaging sealed. Keep the product out of reach of children.
Observe the use by date on the packaging.

9. Advice on Cleaning and Disposal
The test for vaginal fungal infections is only intended for one single use. The test can be disposed of with the normal household waste. Separate the individual parts of the packaging according to material and dispose of them in the correct recycling system.
Observe the local and regional guidelines for disposal.

10. Material Composition
The test cassette contains dried antibodies and latex particles applied on the test strip, and ready-to-use extraction/running buffer sealed in the upper part of the cassette.

11. Explanation of symbols used

	Manufacturer
	Attention
	Observe Instructions for Use
REF	Article Number
LOT	Batch designation
	Store between 2 - 30° C. Do not freeze.
STERILE R	Sterilised by irradiation
	Only for single use.
	Content sufficient for 1 test
IVD	Only for in vitro diagnostics (For external use only)
	To be used by: (Year and Month)

Questions and answers
What shall I do if the test result was negative but the complaints persist?
Please consult a physician. He/she could diagnose other infectious diseases.

What should I do if the test result is positive?
Please see a doctor. He or she will confirm the diagnosis and discuss the further treatment with you.

If the test result is positive: must both lines be equally strong?
Not necessarily. It is irrelevant which line is stronger. Even a weak second line indicates a positive test result.

How important is it to observe the window of time in reading the result?
It is important that you read the result after 10 minutes and not longer than 20 minutes.

How do I know if I used the smear swab correctly?
You will know as soon as the control line becomes visible. The test is sufficiently sensitive to trace even small amounts of vaginal yeast.

The test kit fell down. Can I still use it?

You can use the test as long as all components are intact. If one of the test components is damaged, however, the test can no longer be applied.

Is there a preferable time of day for applying the test?
No. The test can be performed at any time of day.

What shall I do if parts of the liquid run out of the sealing cap when I remove the aluminium cover or afterwards?

The test is invalid in these cases. The amount of liquid in the sealing cap is essential for the correct function of the test.

Can I still use the test if the unpacked smear swab falls down without the protective foil?

No. Due to hygienic reasons we advise not inserting a not-intact smear swab into the vagina.

I turned the sealing cap incorrectly before I removed the aluminium cover or before I put the smear swab into the sealing cap:

Is it possible to bring it back into the basic position and to perform the test?

No. As soon as the upper part of the test has been turned you cannot change it back. The test is now invalid.

Specialised information on how the test works

Vulvovaginal candidiasis (VVC) is thought to be one of the most common causes of vaginal symptoms. Approximately 75 % of women will be diagnosed with a yeast infection at least once in their lifetime. Recurrent infections are experienced by 40 - 50 % of women, and 5 % are estimated to develop chronic yeast infections.

Yeast infections are more often misdiagnosed than other vaginal infections (trichomonas infection and bacterial vaginosis). Symptoms of VVC may include acute itching, vaginal pain, inflammation, a rash on the outer labia of the vagina and genital burning that may increase during urination. However, these symptoms are non-specific. Clinicians should keep in mind that a broad variety of infectious and non-infectious diseases can cause a similar array of symptoms. To obtain an accurate diagnosis, a thorough examination is necessary.

The following standard tests should be carried out in women complaining of vaginal symptoms: measurement of the vaginal pH, saline and 10 % potassium hydroxide microscopy.

Microscopy is the mainstay in the diagnosis of VVC. Recent studies show that, in academic settings, microscopy has a sensitivity of 50 % at best, and thus will miss a substantial percentage of women with symptomatic VVC. To increase the accuracy of diagnoses, yeast cultures have been advocated by some experts as an adjunctive diagnostic test. These cultures are associated with cost and inadequate, and they have the disadvantage that it may take up to one week to get a positive result. Inaccurate diagnoses of yeast infections may delay treatment and cause more serious vaginal diseases.

The test for vaginal yeast infections

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