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Genauigkeit: 90,8 %

### Literatur:

Th. Heinze et al., Determination of vaginal pH by pH indicator strip and by pH micro electrode, Journal of Perinatal Medicine - Official Journal of the WAPM, Band 17, Heft 6, Seiten 477–479

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
## **GB** Instructions for use ELANEE vaginal pH self test

Dear Customer,

Thank you for choosing our ELANEE vaginal pH self test! These instructions should help you to get the best use out of the product. Please read the instructions carefully – particularly the safety and application instructions – before using this product for the first time.

We will be glad to answer any questions and to receive any suggestions you might have.

*Your ELANEE Team*

 Please observe the instructions for use and the packaging and keep them in a safe place.

### **1. Applicability**

These instructions for use apply to the following products:

**REF** 722-00 ELANEE vaginal pH self test, 2 pcs    **REF** 723-00 ELANEE vaginal pH self test, 20 pcs

### **2. Contents**

<b>ELANEE vaginal pH self test, 2 pcs</b>	2 test applicators, individually packaged in plastic wrappers with printed colour scale, 1 instruction leaflet
<b>ELANEE vaginal pH self test, 20 pcs</b>	20 test applicators, individually packaged in plastic wrappers with printed colour scale, 1 instruction leaflet

The ELANEE vaginal pH self test for home use consists of a vaginal applicator about 8 cm long with an integrated pH test strip (MERCK pH indicator strip, pH 4.0–7.0)

### **3. Intended purpose and appropriate use**

The ELANEE vaginal pH self test is a fast solution for testing vaginal pH values. As an in vitro diagnostic agent, this over-the-counter device is suitable for determining bacterial infections in vaginal secretion.

The test result may indicate a shift in vaginal equilibrium. If a shift is detected in good time, corrective measures may be taken at an early stage. Disruptions to the vaginal flora promote bacterial infections. The interpretation of the test result (**bacterial infection, fungal infection or Trichomonads**) in combination with your own symptoms is to be found on pages 28 and 29.

**As soon as the test detects initial indications of a divergence from normal values, a full medical checkup is recommended.**

The ELANEE vaginal pH self test is not suitable for any purposes other than the listed intended ones.

### **4. Field of application**

The indicator paper requires moistening with vaginal secretion. The device is therefore intended for female users.

The ELANEE vaginal pH self test is designed for use in the following cases:

- Bacterial vaginosis
- After anti-infective treatment (e.g. with antibiotics)

- As part of precautionary checkups
- During pregnancy
- When there is a known risk of premature birth
- After sexual contact with a new partner
- For complaints in the genital area (burning, itching, unpleasant odour, discharge)
- Before and after surgery in the vagina

## **5. How it works**

The ELANEE vaginal pH self test is a reliable, pain free method for testing the vaginal pH value. The applicator's pH strip changes colour immediately when it comes into contact with vaginal secretion. This colour is the test result corresponding to the matching pH value on the scale. The regular use of the ELANEE vaginal pH self test for determining vaginal pH values serves as a reliable early warning of disruption to the vaginal flora.

## **6. Medical background and restrictions in use**

The vagina's protective system can only function properly when the vaginal pH value is acidic (between 3.8 and 4.4). This serves effectively to combat incipient and pathogenic infections of the vagina. There are many factors involved in the regulation and function of the vaginal microbial ecosystem and the maintenance of its natural protective role. One of the principal factors is lactobacilli. A healthy vaginal flora is therefore the most important and most natural protection against vaginal diseases.

Certain conditions may temporarily alter the vaginal pH value and so falsify the results of this test. Please note the following restrictions when testing and interpreting the results.

When to test:

- at least 12 hours after sexual intercourse without a condom (sperm may distort the test result)
- at least 12 hours after the application of vaginal medication (pessaries, creams, gels, etc.)
- not during and up to 3 to 4 days following the end of menstruation
- at least 15 minutes after urination (traces of urine may falsify the test result)
- not after washing or showering the external genitals.

Under certain circumstances, the vaginal pH value may be affected by the following factors:

- Bacterial vaginosis
- Mixed bacterial infections
- Sexually transmitted diseases
- Waters breaking early (P-PROM)
- Oestrogen deficiency
- Postoperative wound infections
- Excessive intimate care
- Treatment with antibiotics

### False positive:

A false positive is a pH value result greater than that of the tested vaginal secretion. This can be put down to contamination e.g. with urine or amniotic fluid.

### False negative:

A false negative is a pH value result less than that of the tested vaginal secretion. This can be put down to residual products, e.g. creams, gels, vaginal suppository, that briefly acidify the vaginal secretion and are still in the vagina shortly before the test.


## **7. Usage and interpretation**

Outside of pregnancy, recommended use is once a month or when infection is suspected.

During pregnancy, the test should be applied two to three times a week. This verifies that the vaginal pH value is within limits, e.g. between visits to the gynaecologist. Bacterial vaginosis is diagnosed in about 15 to 20 % of pregnant woman. By regularly testing their vaginal pH value, pregnant woman can reduce their risk of infection.

The procedure must be followed exactly for reliable results. Testing the vaginal pH value is easy and pain free.


### Usage

1. Wash your hands before applying the test.
2. Open the box, take out a test device, and open the plastic wrapper. Do not remove the applicator until just before the test.
3. Hold the applicator by its grip zone between thumb and forefinger (see Fig 1).  
The thumb should be on the side of the applicator without the pH measuring zone.
4. Avoid contact between the pH measuring zone and other fluids.
5. Sit down with your legs spread, or stand with one foot on the toilet seat or other stable surface.
6. With your free hand, carefully open the outer lips (outside of the vaginal opening).


7. Carefully insert the applicator like a tampon 1 to 3 cm into the vagina (Fig 2), but no further than the grip ring. Inserting the vaginal applicator is completely pain free.
  8. Carefully press the inserted applicator against the rear vaginal wall (towards the rectum) for about 10 seconds so that the pH test area is well moistened.
- ⚠ CAUTION! Adequate moistening is crucial for the test results.**
9. Now carefully remove the applicator from the vagina.

10. Read the result immediately after removing the applicator from the vagina when the test area is still moist (Fig 3). Compare the colour of the pH measuring area against the colour scale, and identify the colour closest to your result.
11. Note the result in the test diary. (See table in the instructions.)
12. After use, dispose of the applicator in household waste.

## Interpreting the test results

or

**The pH value of your vagina is normal (3.8 to 4.4).**


Vaginal pH values between 3.8 and 4.4 are in the normal/healthy range. The indicator pad turns a yellowish colour.

If the value is less than 4.4, but there are other indications of infection, seek medical advice.

A **fungal infection** is probable when you notice the following additional symptoms in your genital area:

- Itching
- Burning
- Unusual, whitish, lumpy discharge

**Then seek medical advice.**

The pH value of your vagina is raised (greater than 4.4).

Vaginal values greater than 4.4 are outside of the normal/healthy range. The indicator turns a darker colour, tending towards brown-green-blue.

A **bacterial infection (vaginosis)** is probable when you notice the following additional symptoms in your genital area:

- Itching
- Unpleasant, fishy odour
- Unusual whitish grey, liquid discharge

**Trichomoniasis** is probable when you notice the following additional symptoms in your genital area:

- Pain when urinating
- Unpleasant, musty odour
- Unusual, foam-like yellowish to greenish discharge

If your tested pH value is greatly increased (6.1 or higher), this could indicate **breaking waters** during pregnancy. Consult your gynaecologist or midwife as quickly as possible.

**In the event of raised pH values, do not take any therapeutic measures without first consulting your doctor!**

## 8. Safety instructions

- ⚠ For use exclusively as an in vitro diagnostic agent.
- ⚠ Use each test once only.
- ⚠ Use the test before its expiry date (see date on the box or foil bag).
- ⚠ Do not open the test until just before using.
- ⚠ Do not touch the pH measuring area on the vaginal applicator.
- ⚠ Substances causing an abnormal discolouration of vaginal secretion (menstrual haemoglobin) may affect the legibility of the pH result.
- ⚠ Keep this product out of the reach of children.
- ⚠ The applicator is designed for testing vaginal pH values only. It does not confirm infections.
- ⚠ An acidic pH value is not always adequate in preventing infection. If you detect disease, although the pH value is normal, seek medical advice.
- ⚠ Failing to use the applicator properly may tear the hymen.

## 9. Safekeeping, service life, storage and transport

Please keep the ELANEE vaginal pH self test in the sealed plastic wrapper in a dry place at a temperature between 10 °C and 30 °C. The test should be kept out of the reach of children. If the foil bag is sealed, the product can be used up to the expiry date shown on the packaging or foil bag.

## 10. Cleaning and disposal instructions

The applicator is intended to be used once only. Used vaginal applicators should be wrapped in a tissue when placed in domestic waste. This serves to prevent their misuse. Please separate the individual packaging components and dispose of them in the relevant recycling system. Please be aware of local disposal guidelines.

## 11. Material composition

Each vaginal applicator is made up of the following parts:

<b>Plastic body for inserting in the vagina</b>	ABS plastic
<b>pH measuring area</b>	Test strips of polyethylene terephthalate, and filter paper with specific pH reagents

The strips for the ELANEE vaginal pH self test are nonbleeding, i.e. they do not release colourants or other substances. They may be inserted in the vagina. Do not under any circumstances use pH test strips that release substances (e.g. water test strips).

## 12. Performance features

### Testing range and accuracy:

The testing range of the ELANEE vaginal pH self test is 4.0 to 7.0. Colours range from a soft orange, to yellow and green, to blue for a pH range of 4.0 to 7.0. The colour closest to the colour on the pH scale represents the pH value. Values less than 4.0 are shown as 4.0. Values higher than 7.0 are shown as 7.0. The test is therefore sufficiently accurate for differing between samples with normal and those with atypical pH values.

### Diagnostic sensitivity and specificity:

The suitability of the pH indicator for testing vaginal pH values was verified in the publication by Heinz et al. (1989) which compared the pH results with the values returned by a microelectrode placed directly in the vagina. Samples of up to pH 4.4 were classified as positive samples and described as the normal vaginal environment. Samples of greater than pH 4.4 were classified as negative. These samples describe an atypical vaginal environment with an increased risk of infection.

Sensitivity: > 99.9 %

Specificity: 87 %

Positive predictive value (PPV): 75.6 %

Negative predictive value (NPV): > 99.9 %

Accuracy: 90.8 %

### Literature:

Th. Heinze et al., Determination of vaginal pH by pH indicator strip and by pH micro electrode, Journal of Perinatal Medicine - Official Journal of the WAPM, Volume 17, No. 6, pp 477–479

### Interference test:


Acetaminophen, acetic acid, ampicillin, ascorbic acid, atropine, beta-Hydroxybutyric acid, bilirubin, caffeine, creatinine, EDTA, ephedrine, estradiol, estriol, ethanol, glucose, glycerol, haemoglobin, human serum albumin, Ibuprofen, lactic acid, phenethylamine, pregnanediol, progesterone, riboflavin, salicylic acid, tetracyclin, urea, and uric acid were tested for their interference (superimposition) with results. With the exception of haemoglobin, none of these substances falsify the results for pH values 4.0, 4.4, 4.7, 5.0, 5.5, 5.8, and 7.0.

### Limits

The effects of medication and other metabolic products on the ELANEE vaginal pH self test are not known in all cases. In cases of doubt, the test should be applied again after medical consultation.

It should be borne in mind that there are certain conditions that can temporarily affect the pH value (see restrictions in use and interference), and that the user must make sure that these conditions do not arise.

### 13. Explanation of the symbols used

	Manufacturer		Not intended for repeated use
	Warning		Contents adequate for 2 tests (2 pcs)
	Follow instructions for use		Contents adequate for 20 tests (20 pcs)
	Item number		To be used by (year month)
	Batch designation		Protect from moisture
	Store at 10 °C - 30 °C		In vitro diagnostic agent
	Keep out of direct sunlight		

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