ELANEE®

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---

### Instructions for Use

#### ELANEE seat cushion for pelvic floor training

##### Dear customer,

Thank you for choosing our ELANEE seat cushion! These instructions will help you utilise our product as well as possible. Please read the instructions for use thoroughly before your first training.

We are available if you have any questions or suggestions and wish you much success with your pelvic floor training!

Your *ELANEE Team*

Please observe and keep the instructions for use and the packaging!

#### 1. Field of application

These instructions for use apply to the following product:

708-00 ELANEE seat cushion for pelvic floor training

#### 2. Product content

- 1 seat cushion for pelvic floor training
- 1 hand pump
- 1 instructions for use
- 1 Exercise program

#### 3. Purpose and intended use

The pelvic floor is an approximately palm-sized, flat and elastic group of muscles between pubic bone, tailbone and the ischial tuberosity. The pelvic floor muscles consist of three muscles layers. They perform the closing function of the urethra and anus, support the internal organs in the abdomen and affect posture. The pelvic floor is crossed by the urethra, vagina and anus.

Due to the hormone changes during pregnancy and childbirth, the climacteric period or menopause, due to severe physical strain or overweight, overexpansion or weakening of the pelvic floor can occur. The consequences can be bladder weakness or a malposition of the organs in the abdomen.

The pelvic floor can however be trained just as every other muscle by intense training. This way, you can correct these discomforts or avoid them preventively. Pelvic floor training is particularly suited:

- During pregnancy and after childbirth for preventing pelvic floor weakness and for supporting the regeneration
- For weakening, correcting and preventing bladder weakness (due to stress incontinence) with women (particularly in the climacteric period and menopause)
- For improvement, if the uterus (Descensus uteri) or other pelvic organs have descended

#### **4. Functionality**

With the ELANEE seat cushion for pelvic floor training you can train your pelvic floor muscles on the side in everyday life. The operating principle of the seat cushion is “active sitting”. To do that, the seat cushion is simply placed on any seating furniture. Unpadded surfaces are particularly suited.

The cushion is filled with air and creates a slight instability in the pelvis during sitting. The pelvic floor muscles counterbalance this instability and are therefore constantly strained and stimulated during sitting. The constant movement trains the pelvic floor muscles. The ripped side of the cushion stimulates blood circulation in the buttocks at the same time.

The ELANEE seat cushion for pelvic floor training is a perfect addition to pelvic floor gymnastics or training with the pelvic floor training aids (e.g. ELANEE pelvic floor training aids phase I or phase II).

#### **Other positive effects of “active sitting” are:**

- Training of weakened pelvic floor muscles and preservation of existing muscle strength
- Positive effect on upright and healthy posture
- Training of the back muscles

- Helps in prevention and treatment of back problems
- Helps in preventing disc prolapse

## 5. Application

1. Check if the cushion is in perfect condition.
2. Place the seat cushion with the ripped side upwards on solid and level seating.
3. Sit on the cushion.
4. Sit on the seat cushion as long as it is comfortable for you and / or do 10 minutes of exercises on it every day.
5. If the cushion loses air after using it, simply pump it back up with the hand pump provided. To do that, screw the needle onto the pump and then insert it into the intended valve of the cushion. Now press the pump until the desired air volume is reached. In order to remove air from the cushion, simply push the needle into the valve without the pump and gently squeeze the cushion until you have the volume of air you want.

**Tip:** In the provided exercise programme you will find exercise suggestions for effective pelvic floor training with the seat cushion.

## 6. Safety information

Only use the seat cushion on solid and level seating!

Avoid contact with sharp or pointed objects in order to prevent damaging the cushion!

Do not use abrasive cleaning agents or disinfectants containing alcohol or chlorine and do not use hot air for drying!

## 7. Cleaning and care instructions

- The seat cushion can be cleaned with a mild cleaning agent and a moist cloth.
- Please separate the individual components of the packaging according to different materials and recycle them appropriately. Please observe the local and regional disposal regulations in the process.

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