ELANEE®

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Achtung!

Trocken lagern

## Instructions for use ELANEE Pilates Ring

### Dear Customer,

Thank you for purchasing our ELANEE Pilates ring! These instructions for use shall help you to use the product optimally.

Please read the instructions for use carefully- in particular the „safety advice” - before using this product for the first time.

We will be glad to answer any questions and receive any suggestions you might have.

Your ELANEE-Team

Read and keep instructions for use and packaging!

### 1. Area of Application

These instructions for use apply to the following product:

**REF** 709-V1 ELANEE Pilates ring

### 2. Contents

- 1 ELANEE Pilates ring
- 1 set of instructions for use
- 1 exercise programme

### 3. Purpose and Intended Use

The ELANEE Pilates ring is a simple and effective training device for a large number of exercises. The ring can be used for active full body training to strengthen the pelvic floor, stomach and back muscles. This helps to stabilise the spine and pelvic floor, while training co-ordination and balance.

The ring is also suitable for targeted strength and muscle-building training in the chest, shoulder, arm and leg areas.

The training device is for people of any age and can be used by anyone, from beginners to advanced practitioners. It can be used for individual or group training.

The ELANEE Pilates ring has been conceived exclusively for use in dry indoor areas.

The ring should only be used for therapeutic purposes after consultation, and with specialist support (e.g. with a physiotherapist).

The ELANEE Pilates ring is not suitable for any purposes other than the listed intended ones.

#### 4. Area of Application

The ring is an effective training device and can be used by men and women during Pilates or full body training - regardless of whether the user is a beginner, or an advanced practitioner.

The ELANEE Pilates ring is especially suited to women with weakened pelvic floor, stomach and back muscles.

#### 5. Functionality

The ELANEE Pilates ring is placed between the arms or legs, depending on the desired muscle group. The ring has limited elastic properties, and provides a marked resistance to pressure from the arms or legs.

The two holders on the ring help with getting the ring into the correct position between the arms or legs. The curved grips provide a better hold and good padding.

The muscles are worked and strengthened by holding and continually pressing on the ring, or alternatively pressing and releasing the ring. Weakened muscle groups such as the pelvic floor, stomach and back muscles can be built up and stabilised, depending on the training position, the intensity of the pressure and the number of repetitions.

#### 6. Application

1. Wear comfortable clothing when training and make sure you have enough space around you. Have a padded sport mat (Pilates mat) ready.
2. Take your time when training
3. Follow the instructions in the attached exercise programme (D, GB, F only).

**Tip:** The length and intensity of the training session can be adjusted to suit your personal level of fitness and comfort.

#### 7. Safety Advice

Please check your ELANEE Pilates ring before each use. Only use a Pilates ring in perfect condition (with no damage or sharp edges).

Pilates ring is under tension!

In order to avoid the possible risk of injury when pressing the ring, the device should not be stored near unsupervised children.

## 8. Storage and Shelf Life

Please store the ELANEE Pilates ring in a dry area, out of sunlight – ideally, in its packaging.

Keep the device out of the reach of children.

The Pilates ring is intended for repeated and long-term use.

## 9. Advice on Cleaning and Disposal

The Pilates ring can be cleaned with a damp cloth after each training session.

The Pilates ring can be disposed of with the usual household rubbish.

Please separate and sort the individual packaging components and put them in the correct recycling systems.

Please be aware of local and regional disposal guidelines.

## 10. Material Composition

The ELANEE Pilates ring consists of:

- Polypropylene (PP)
- Thermoplastic rubber (TPR)
- Nitrile butadiene rubber (NBR)
- Fibreglass

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