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**(GB)** Exercise program with ELANEE Pilates Ring, for exercising the pelvic floor, core and back muscles

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**(GB) Before you begin ...**

The following exercises are particularly helpful if you are aware of your pelvic floor and can actively tense and relax it and if you breathe correctly during the exercises.

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**(GB) Activating the pelvic floor**

Sit directly onto a chair and feel where your sitting bones are. The pelvic floor muscles (and the vaginal and anal openings) are located between these two bones. Activate the muscles by imagining that you are going to suck in these two openings into your body. Making yourself aware of the pelvic floor muscles for the first time is often tricky. It is often useful to take part in a specially-designed course or to use the ELANEE Pelvic Floor Training Aid Phase I.

**Breathing correctly**

We frequently breathe too „weakly“, just in the upper chest area. It is therefore most effective to breathe in through the nose and into the stomach (your ribcage will expand outwards), and then out through the mouth, as if you are trying to blow out a candle.

**Please note**

If you have any physical injuries or health complaints, you should consult your doctor or therapist about whether the following exercises are suitable for you.

**You will need**

To begin, you will need a padded mat, comfortable clothing, the ELANEE Pilates Ring, and of course some time (and motivation).

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## 1 One hundred - Cent

**Starting position - Position de départ**

**(GB)** Lie flat on your back, with your arms lying parallel to the body, and hold the Pilates Ring between your bent legs (hip-width apart), just above the knee.

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
**(GB)** Inhale and raise your head at the same time (look in the direction of your knees). Lift your shoulders and outstretched arms off the floor. Inhale and exhale evenly for 5 seconds each, „pumping“ your arms continuously (quickly move the outstretched arms up and down). Count until you get to 100 pumping movements. On every exhalation, the abdominal wall should flatten and the waist should become narrower.

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**Attention:** Tense your stomach, stretch out your arms, hold your shoulders low.  
**Attention :** Tenir la tension abdominale, bras tendus longs, épaules basses.

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**Starting position - Position de départ**

**(GB)** In a seated position: Place the Pilates Ring between your thighs. Place your feet hip-width apart and parallel to one another. Place your hands under the back of your knees. Point your elbows outwards.

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
**(GB)** Squeeze the ring gently. Breathe in, arch your lower back, bring your chin down to your chest and roll back until your arms are stretched. Exhale, activate your core and pelvic floor muscles and roll up slowly.

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## 3 Corkscrew - Tire-bouchon

**Starting position - Position de départ**

**(GB)** Lie flat on your back, holding the Pilates Ring between your ankles and gently squeeze the ring. Stretch your legs towards the ceiling. Place your hands under the base of the spine or parallel next to the body. Place your fingers underneath your buttocks.

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
**(GB)** Activate your core and pelvic floor muscles. Inhale and move your outstretched legs in the shape of a circle: first towards your nose, then to the right, to the front, and to the left. Meanwhile, control your breathing.

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**Attention:** Press your shoulders/upper arms into the mat. Gently bend your legs if stretching them out causes you problems. Do not circle your legs wider than the mat.

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## 4 Bridge - Pont

**Starting position - Position de départ**

**(GB)** Lie flat on your back, with your arms lying parallel to the body, and hold the Pilates Ring between your bent legs (hip-width apart) just above the knee. Squeeze gently.

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**(GB)** Activate your core and pelvic floor muscles and tense. Exhale and roll up your spine, vertebra by vertebra, from your pelvis to your shoulders. Inhale here and on the exhalation, roll down your spine vertebra by vertebra, until you are lying flat again. Repeat.

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**Attention:** Make sure that your body is straight when lifted. Look up to the ceiling. Press your shoulders/upper arms into the mat. Make sure that your thighs are parallel and knee facing forwards.

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## 5 Recovery position - Position latérale

**Starting position - Position de départ**

**(GB)** In the recovery position: Lie on your side, resting your head on your outstretched arm. Rest your other hand on your ribcage. Your legs should be stretched out parallel to one another. Your shoulders and hips should be aligned. Gently lift your waist. Look directly ahead. Place the Pilates Ring between your lower legs, just above the knee.

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
**(GB)** Activate your pelvic floor and core muscles. Exhale while using the upper leg to press the Pilates Ring into the lower leg. Hold this pressure for three seconds and then relax. Repeat.

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**Attention:** Keep your buttocks completely relaxed. Your hips and shoulders should stay straight.  
**Attention :** Garder le corps entièrement au repos. Les hanches et les épaules restent droits.

## 6 Crab - Crabe

**Starting position - Position de départ**

**(GB)** In the recovery position: Lie on your side, resting your head on your outstretched arm. Rest your other hand on your ribcage. Gently bend your legs. Your shoulders, hips and knees should be aligned. Place the Pilates Ring between your thighs and squeeze gently. Look directly ahead.

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
**(GB)** Hold your legs just above the floor (so that they are not in contact with the mat at all). Activate your deep core and pelvic floor muscles. Exhale and gently squeeze the Pilates Ring. Hold this position for three seconds. Tense your muscles more with each exhalation.

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## 7 Lifting your knees - Lever les genoux

**Starting position - Position de départ**

**(GB)** On all fours: Make sure that your hands are underneath your shoulders and your knees underneath your hips. Open your knees to hip-width apart, keeping your neck long, with your feet up on your toes. Place the Pilates Ring between your thighs and squeeze gently.

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**(GB)** Activate your core and pelvic floor muscles. Hold your back straight, leave your neck long and relax your shoulders. Exhale and gently lift your knees up from the floor. Hold this position for five seconds. Release your knees down to the floor. Repeat.

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**Starting position - Position de départ**

**(GB)** Standing at a wall: Make sure that your head, back and buttocks are touching the wall. Bend your legs until your thighs are parallel to the floor (in a seated position). Open your thighs to roughly hip-width apart, positioning your feet parallel to one another, knees pointing forward, and making sure that your knees do not come over your toes. Place the Pilates Ring between your thighs.

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**(GB)** Activate your core and pelvic floor muscles. Hold your back straight, leave your neck long and relax your shoulders. Look directly ahead. Exhale and gently squeeze the ring. Hold the pressure for 5 seconds.

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**Attention:** Make sure that your back is flat against the wall