ELANEE®

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Artikelnummer

Achtung!

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GB Pilates-Ring – Instructions for use

Dear Customer,

Thank you for choosing our ELANEE Pilates ring!

The Instructions for Use should help you to get the best use out of your product. Please read the instruction manual carefully – particularly the usage and safety instructions – before using this product for the first time.

We will be glad to answer any questions and to receive any suggestions you might have!

Your *ELANEE-Team*

Please follow the instructions and retain the instruction sheet and the packaging!

1. Field of application

These instructions for use are applicable for the following product:

REF 709-V1 ELANEE Pilates ring

2. Product content

- 1 ELANEE Pilates ring
- 1 instructions for use
- 1 exercise programme

3. Purpose and Intended Use

The ELANEE Pilates ring is a simple and effective training device for active full body training to strengthen pelvic floor, stomach and back muscles. The training device is for people of any age and can be used by anyone, from beginners to advanced practitioners. The ELANEE Pilates ring is especially suited to women with weakened pelvic floor, stomach and back muscles.

The Pilates ring is designed for use at home. The ring should only be used for therapeutic purposes after consultation, and with specialist support (e.g. with a physiotherapist). The training device is designed exclusively for use in dry indoor areas. The ELANEE Pilates ring is not suitable for any purposes other than the listed intended ones.

4. Area of Application

The Pilates ring is suitable for targeted strength and muscle-building training in the pelvic floor, stomach and back muscles as well as the chest, shoulder, arm and leg muscles. This helps to stabilise the spine and pelvic floor, while training co-ordination and balance.

5. Functionality

The ELANEE Pilates ring is placed between the arms or legs, depending on the desired muscle group. The ring has limited elastic properties, and provides a marked resistance to pressure from the arms or legs.

The two holders on the ring help with getting the ring into the correct position between the arms or legs. The curved grips provide a better hold and good padding.

The muscles are worked and strengthened by holding and continually pressing on the ring, or alternatively pressing and releasing the ring. Weakened muscle groups such as the pelvic floor, stomach and back muscles can be built up and

stabilised, depending on the training position, the intensity of the pressure and the number of repetitions.

6. Usage

1. Practice in a warm, quiet place where there are no distractions. Make sure there is sufficient space for training and that there are no objects protruding into this space. Have a padded sport mat ready.
2. Wear comfortable clothing when training. Carry out these exercises preferably bare-foot or wearing sports socks. Avoid wearing jewellery or a watch when training.
3. Take your time when training. Before you start the exercises, you should warm up by doing warm-up exercises for approximately 10 minutes.
4. After warming up, follow the instructions in the attached exercise programme (only in German, English and French).
5. The length and intensity of the training session can be adjusted to suit your personal level of fitness and comfort.
6. Finish the training with a relaxation phase of approximately 5 minutes.

Tip: It can be helpful to carry out the exercises in front of a mirror, so that you are able to check your posture more easily.

7. Safety Instructions

-  Please check your ELANEE Pilates ring before each use. Only use a Pilates ring in perfect condition (free from cracks, damage or sharp edges).
-  Pilates ring is under tension! Do not overload it, otherwise there is a risk of injury.
-  Do not use the Pilates ring for any other purpose. Do not hang from the Pilates ring, do not sit or stand on it and do not use it to hang up other objects.
-  When using during pregnancy, shortly after birth, if you have any existing injuries, after recent operations or if there are other impairments, such as, for example, inflammatory illnesses of the tendons or joints, you have a pacemaker, orthopaedic conditions, consult with your doctor to agree a training programme using the device, if necessary.
-  Stop the training immediately and contact your doctor if any of the following symptoms should occur: chest pains, severe shortness of breath, dizziness, nausea, pains in the joints and muscles.

 The Pilates ring is not a toy. The training device must be kept out of reach of children, in order to avoid any possible risk of injury.

 Do not train if you are distracted, as this could mean there is a risk of injury.

8. Safekeeping, service life, storage and transport

Please store the ELANEE Pilates ring in a dry area, out of sunlight – ideally, in its packaging.

Keep the device out of the reach of children.

The Pilates ring is intended for repeated and long-term use.

9. Cleaning and disposal instructions

The Pilates ring can be cleaned with a damp cloth after each training session.

The Pilates ring can be disposed of with the usual household rubbish.

Please separate and sort the individual packaging components and put them in the correct recycling systems.

Please be aware of local and regional disposal guidelines.

10. Material Composition

The ELANEE Pilates ring consists of:

- Polypropylene (PP)
- Thermoplastic rubber (TPR)
- Nitrile butadiene rubber (NBR)
- Fibreglass

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