ELANEE®

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## ELANEE Beckenboden-Trainingshilfe

### Phase II

|           |   |    |
|-----------|---|----|
| <b>D</b>  | Beckenboden-Trainingshilfe Phase II –<br>Gebrauchsanweisung .....                                   | 2  |
| <b>GB</b> | Pelvic Base Training Aid Phase II – Instructions for use .....                                      | 7  |
| <b>F</b>  | Accessoire de rééducation du plancher pelvien niveau II –<br>Notice d'utilisation .....             | 12 |
| <b>I</b>  | Ausilio per l'allenamento dei muscoli del diaframma pelvico<br>fase II – Istruzioni per l'uso ..... | 18 |
| <b>NL</b> | Bekkenbodem Trainingshulp Fase II – Gebruiksaanwijzing...   | 24 |
| <b>ES</b> | Ayuda de entrenamiento del suelo pélvico fase II<br>Instrucciones de uso. ....                      | 29 |
| <b>RU</b> | Кегель для укрепления мышц тазового дна фаза II –<br>Инструкция по применению.....                  | 35 |

**LOT** siehe Verpackung / see packaging / voir l'emballage/  
vedi confezione / zie verpakking / véase el embalaje / См. упаковку

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Artikelnummer

---

## **(GB) Pelvic Base Training Aid Phase II – Instructions for use**

**Dear customer,**

Thank you for choosing our ELANEE Pelvic Base Training Aid Phase II! The Instructions for Use should help you to get the best use out of your product. Please read the instruction manual carefully – particularly the usage and safety instructions – before using this product for the first time. We will be glad to answer any questions and to receive any suggestions you might have!

*Your ELANEE team*

Please follow the instructions and retain the instruction sheet and the packaging!

Made in Germany

### **1. Field of application**

These instructions for use apply to the following product:

705-00 ELANEE Pelvic Base Training Aid Phase II

### **2. Product content**

- 1 Training aid, 1 Instructions for use

### 3. Purpose and intended use

The ELANEE Pelvic Base Training Aid Phase II is a product for vaginal use. It is used for long-term follow-up training of the pelvic floor muscles when the pelvic floor can already be perceived and basic pelvic floor strength has been achieved (e.g. with the ELANEE Pelvic Floor Training Aid Phase I).

The ELANEE Pelvic Base Training Aid Phase II is not suitable for any purposes other than the listed intended ones.

### 4. Areas of application

- During a **normal pregnancy**, you can use the Pelvic Base Training Aid Phase II as long as you feel comfortable.
- **After giving birth**, you can start exercising as soon as the cervix has closed (approx. six weeks postpartum) and your pelvic floor muscles have regained some strength. If in doubt, please ask your midwife or doctor.
- You can also use the Pelvic Base Training Aid Phase II to train your pelvic floor **during the climacteric period**.

### 5. Functionality

The ELANEE Pelvic Base Training Aid Phase II is composed of two tubular bodies that are interconnected by a wide strip. Both tubular bodies contain a ball floating freely. The internal balls swing against the outside wall of the tubular body with every motion and generate slight vibration. This vibration is directly transferred to the pelvic base muscles. This stimulates and activates the muscles.

The large supporting surfaces of the lower part help transfer the vibrations directly to the pelvic base muscles. The 14 imprinted stimulation islands facilitate additional stimulus to the muscles.

The wide strip between both halves enhances the swinging of both balls and the training effect is intensified by this.

The ELANEE Pelvic Base Training Aid Phase II is fully inserted into the vagina. For easier insertion, the upper part is smaller than the lower part. The training aid is completely silent and cannot be seen from the outside. The muscles are trained simply by wearing the training aid, making it easy to incorporate into normal daily routines. The pelvic floor muscles are stimulated with every movement.

A strong pelvic floor can also have a positive effect on sexual sensation.

## **6. Usage**

1. Wash the Pelvic Base Training Aid Phase II thoroughly before the first use and make sure that conditions are clean (e.g. wash hands).
2. Insert the small upper part first. The use of medical lubricating gel can be helpful when training with the ELANEE Pelvic Base Training Aid Phase II, in particular if your vagina is dry..
3. Insert the training aid slowly, completely and approximately as far as you would insert a tampon. Both crescent-shaped supporting surfaces should be right and left during insertion. This way, they rest directly on the pelvic base muscles.
4. Try in the beginning how long you feel wearing them is comfortable and increase the length according to your well-being.
5. Conduct the training daily so that your pelvic base muscles remain active.

## 7. Safety Instructions

-  Thoroughly clean the training aid before first use and after each subsequent use and take care to ensure clean conditions (e.g. wash your hands).
-  Please check your training aid before each use. Only use a product if it is in perfect condition. A training aid with a scratched or damaged surface could cause injuries.
-  The pelvic base training aid should not be used in case of vaginal infections, inflammation, injuries or painful conditions or during menstruation or high-risk pregnancy.
-  When other complaints (such as a dropping of the pelvic organs) exist, you should have other causes than a weakness of the pelvic floor ruled out in advance by a doctor before you start with the training.
-  In cases of a known hypersensitivity and allergies to one of the named ingredients, you should not use the aids or only use them after consultation with a doctor.
-  The pelvic floor training aid is not suitable for use in virginity.
-  The wearing of cervical caps, contraceptive rings, diaphragms or intrauterine devices during training are possible. However, you must check the correct positioning before and after each training session (if in doubt only use after consultation with a doctor).
-  The training aid must not be worn during X-rays, a CT scan (computer tomography) or an MRI scan (magnetic resonance imaging).
-  Do not use steam sterilisation, dry heat sterilisation or

boiling to clean the aid as this can lead to distortion of the material.

## 8. Safekeeping, service life, storage and transport

There are no special instructions for keeping, storing and transporting the training aid.

## 9. Cleaning and disposal instructions

- After each use, clean the training aid with warm water and soap or disinfect it with a commercial skin disinfectant.
- Before using the training aid again, any soap or disinfectant residues should be thoroughly rinsed off with water.
- The ELANEE Pelvic Base Training Aid Phase II can be disposed of with household waste.
- Please separate the individual components of the packaging into groups of identical materials and dispose of them in the appropriate recycling system.
- When doing this please note the local and regional waste disposal regulations.

## 10. Material composition

The ELANEE Pelvic Base Training Aid Phase II consists of:

- glass-fibre reinforced polypropylene coated with a thermoplastic elastomer
- two rubberised steel balls
- polyamide yarn coated with polyurethane

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