ELANEE®

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- Remove a test strip from the plastic tube and reseal the plastic tube.
- Ensure that no discolouring is shown on the unused test strip.
- ⚠️ ATTENTION! Do not touch the test fields on the test strip! If the test pads on the test strip show any discolouring, use a new test strip.
- Dip the test strip briefly (max. 1 sec.) into the urine so that all the test fields are moistened with urine.
- Wipe any excessive urine off on the rim of the cup or wipe it off with a liquid-absorbing paper tissue. When doing this, make sure that you only touch the edges of the test strip lightly with the absorbing paper.
- ⚠️ ATTENTION! Do not touch the test fields!
- Read the test results carefully at 60 seconds (leukozytes 120seconds) in a good light and with the test area held near the appropriate color chart on the pouch label. Changes in color that appear only along the edges of the test pads or after more than 2 minutes have passed are of no diagnostic significance.

Expected values:**Glucose:**

The kidney normally excretes small amounts of glucose. Concentrations of 50mg/dl or above are not considered to be normal and are evaluated as a positive result.

The presence of traceable quantities of glucose in the urine often suggests diabetes mellitus. But other causes such as pregnancy or congenital malfunction of the urinary system may also be responsible for glucose in the urine. However, non-existence of glucose in the urine does not exclude an illness.

Protein:

Normal urine specimens ordinarily contain some protein. A continuously high concentration of more than 20 mg/dl are not considered to be normal and evaluated as a positive result. The presence of a continuously high amount of protein in the urine can be an indication of kidney or urinary tract diseases. The evidence of protein in the urine alone is, however, not an unambiguous diagnosis for an illness. There may be harmless causes for protein being present in the urine, such as physical stress.

Nitrite:

Normally no nitrite is detectable in urine. If the test pad changes to pink-red, this indicates that there is nitrate in the urine and the result is positive. In this way it is indirectly proven that nitrite-forming germs are present in the urine. The most frequent pathogens of urinary infections such as Escherichia coli and most of the other urinary pathogenic germs change nitrate that is absorbed with food into nitrite. Even just a weak change to pink indicates a considerably increased amount of bacteria. The test is thus suitable for early detection and for preventing a relapse of urinary infections.

Leukozyten:

Normally no leukocytes are detectable in urine. Concentrations of 20 leukocytes/ μ L or above are not considered to be normal and evaluated as a positive result. Leukocytes are increasingly excreted in the event of kidney infections and the upper urinary tract.

Blood:

Blood may often be found in the urine of menstruating females. The blood concentration in the urine of healthy persons is very low (0-5 erythrocytes/ μ L). Concentrations of above 5-10 Erythrocytes/ μ L are not considered to be normal and are evaluated as a positive result. Haemoglobin is the red blood pigment in the red blood cells (erythrocytes). Separate colour scales are shown on the plastic tubes for erythrocytes and haemoglobin. Single or closely collected green spots on a yellow test patch indicate intact erythrocytes. Haemoglobin or damaged erythrocytes are indicated by a uniform green colouring of the test pad. The main causes for blood to be present in urine are kidney and urinary tract diseases. But serious physical stress, burns or infections (e.g. bladder infection) may also be responsible.

General information:

First morning urine:
This is the urine, which has been collected overnight in the bladder and is obtained in the morning before the breakfast and other activities.

Second morning urine:

The second morning urine is spontaneous urine collected before midday and after the first morning urine.

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