ELANEE®

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EN **Urine-Control Test⁵ – Instructions for Use**

Dear Customer,

Thank you for purchasing our ELANEE Urine–Control Test⁵.

These instructions for use are to help you to use the product properly.

Please read the instructions through carefully- in particular the “safety advice” and “usage” - before using this product for the first time. We will be glad to answer any questions and receive any suggestions you might have.

Your ELANEE Team

⚠ Please read through and keep the instructions for use and packaging!

1. Field of application

These instructions for use apply to the following product:

REF 726-00 Urine–Control Test⁵

2. Product content

- 10 Test strips
- 1 Plastic tube with colour chart and lid with desiccant
- 1 instruction for use

Additional special equipment:

- Timer (stopwatch)
- Absorbent paper (tissue, gauze)
- Clean dry container (tube, cup)

3. Purpose and intended use

The ELANEE Urine-Control Test⁵ is a quick test, which is easily and simply applied for your own use at home, and which shows various parameters in the urine. Test result may provide information regarding the status of carbohydrate metabolism, kidney function and urinary tract infection. It is measured by comparison of test paper attached to a plastic strip with the color chart printed on the pouch label. The result is available after 1 minute (leukocytes after 2 minutes).

If any symptoms or signs of positive results are obtained after doing this test, consulting a doctor for a complete medical history and physical examination is recommended.

The test is not suitable for any other purposes than the ones named above.

4. Areas of application

Test strips are dip-and-read test strips for In Vitro Diagnostic Use only for testing above items in urine. **Glucose [Sugar], Blood, Protein, Nitrit, Leukozytes (white blood cells).**

5. Functionality

Urine to be tested is absorbed in each test item (portion) on the strips and color reaction is developed of each test item (portion) on the strip by means of chemical & enzyme reactions.

6. Usage

Specimen collection and preparation:

Use only clean dry container to collect urine test it as soon as possible. If testing cannot be done with an hour after voiding, refrigerate the specimen immediately and let it return to room temperature before testing. The improperly stored (stored over 4hours at RT and refrigerator) urine specimens may give inaccurate results.

Visual test procedure:

The procedure must be followed exactly to achieve reliable results.

1. Check the expiry date on the test unit. Do not use after the expiry date. Check whether the colour scale on the plastic tube is undamaged so that you can read the result precisely.
2. Ensure clean conditions (e.g. wash your hands) and have a watch with a second hand available to record the time, as well as a clean container for the urine sample.
3. Use fresh urine in a clean and dry container (e.g. plastic cup)
4. Remove a test strip from the plastic tube and reseal the plastic tube.
5. Ensure that no discolouring is shown on the unused test strip.

⚠ **ATTENTION!** Do not touch the test fields on the test strip! If the test pads on the test strip show any discolouring, use a new test strip.

6. Dip the test strip briefly (max. 1 sec.) into the urine so that all the test fields are moistened with urine.
7. Wipe any excessive urine off on the rim of the cup or wipe it off with a liquid-absorbing paper tissue. When doing this, make sure that you only touch the edges of the test strip lightly with the absorbing paper.

- ⚠ **ATTENTION!** Do not touch the test fields!
8. Read the test results carefully at 60 seconds (leukocytes 120seconds) in a good light and with the test area held near the appropriate color chart on the pouch label. Changes in color that appear only along the edges of the test pads or after more than 2 minutes have passed are of no diagnostic significance.

Expected values:

Glucose:

The kidney normally excretes small amounts of glucose. Concentrations of 50mg/dl or above are not considered to be normal and are evaluated as a positive result.

The presence of traceable quantities of glucose in the urine often suggests diabetes mellitus. But other causes such as pregnancy or congenital malfunction of the urinary system may also be responsible for glucose in the urine. However, non-existence of glucose in the urine does not exclude an illness.

Protein:

Normal urine specimens ordinarily contain some protein. A continuously high concentration of more than 20 mg/dl are not considered to be normal and evaluated as a positive result. The presence of a continuously high amount of protein in the urine can be an indication of kidney or urinary tract diseases. The evidence of protein in the urine alone is, however, not an unambiguous diagnosis for an illness. There may be harmless causes for protein being present in the urine, such as physical stress.

Nitrite:

Normally no nitrite is detectable in urine. If the test pad changes to pink-red, this indicates that there is nitrate in the urine and the result is positive. In this way it is indirectly proven that nitrite-forming germs are present in the urine. The most frequent pathogens of urinary infections such as Escherichia coli and most of the other urinary pathogenic germs change nitrate that is absorbed with food into nitrite. Even just a weak change to pink indicates a considerably increased amount of bacteria. The test is thus suitable for early detection and for preventing a relapse of urinary infections.

Leukocyte:

Normally no leukocytes are detectable in urine. Concentrations of 20 leukocytes/µL or above are not considered to be normal and evaluated as a positive result. Leukocytes are increasingly excreted in the event of kidney infections and the upper urinary tract.

Blood:

Blood may often be found in the urine of menstruating females. The blood concentration in the urine of healthy persons is very low (0-5 erythrocytes/µL). Concentrations of above 5-10 Erythrocytes/µL are not considered to be normal and are evaluated as a positive result. Haemoglobin is the red blood pigment in the red blood cells (erythrocytes). Separate colour scales are shown on the plastic tubes for erythrocytes and haemoglobin. Single or closely collected green spots on a yellow test pad indicate intact erythrocytes. Haemoglobin or damaged erythrocytes are indicated by a uniform green colouring of the test pad. The main causes for blood to be present in urine are kidney and urinary tract diseases. But serious physical stress, burns or infections (e.g. bladder infection) may also be responsible.

General information:

First morning urine:

This is the urine, which has been collected overnight in the bladder and is obtained in the morning before the breakfast and other activities.

Second morning urine:

The second morning urine is spontaneous urine collected before midday and after the first morning urine.

Spontaneous collected urine:

Can be collected at any time. This is the most convenient form but does, however, present the largest deviations.

Midstream urine:

Midstream urine is that part of the urine portion that is taken in the middle of passing urine. The first and last part of the urine cannot be used as samples, as these are contaminated by microorganisms, vaginal secretion, or deposits in the urethra.

False-Positive Results:

A false positive result indicates the presence of substances although there are no such substances in the urine sample.

False-Negative Results:

A false negative result indicates the non-existence of substances although there are substances present in the urine.

Parameter	Recommended urine sample	False-Positive Results	False-Negative Results
Glucose	Midstream urine from the second morning urine (two hours after a carbohydrate meal)	Improper storage, chlorine bleach or other strong oxidizing	High SG [>1.020] with high pH urine and ascorbic acid (more than 40mg/dL)
Protein	Midstream urine from the first morning urine	Strongly basic urine (pH 9) or turbid urine specimens	Other proteins than “albumin” are present
Nitrite	Midstream urine from the first morning urine	Bacterial contamination of the urine specimen, Medication	Ascorbic acid ($>30\text{mg/dL}$) with low level of nitrite containing ($<0.03\text{mg}$), insufficient incubation time in the bladder, insufficient dietary nitrate (e.g. fasting)

Leukocytes	Midstream urine from the first morning urine or spontaneous collected urine	Contamination of the specimen by vaginal discharge	Large amounts of ascorbic acid, high concentration of glucose, high specific gravity, high level of albumin, high concentration of formaldehyde or presence of blood
Blood	Midstream urine from the first morning urine	Microbial peroxidase associated with urinary tract infection. Blood as a contaminant from menstruation	Ascorbic acid concentrations [$>30\text{ mg/dL}$]

7. Safety instructions

- ⚠ For in vitro diagnostic use only. Do not consume.
- ⚠ Each test strip is for a single use.
- ⚠ Do not open pouch until ready to use.
- ⚠ Do not touch test areas of urine reagent strips.
- ⚠ Keep out of reach of children
- ⚠ Do not remove desiccant from box.
- ⚠ The desiccant in the pouch is non-toxic for your health but in case it is inadvertently ingested, flush down it with plenty of water.

8. Safekeeping, service life, storage and transport

Please store the ELANEE Urine–Control Test⁵ in a dry place at room temperature of 2°C to 30°C, keeping the packaging sealed. Do not store the strips in a refrigerator or freezer. Store away from moisture and light.

9. Cleaning and disposal instructions

The ELANEE Urine–Control Test⁵ is only intended for one single use.

Used test strips should be wrapped by tissue (or gauze) and then disposed into waste box to avoid misuse of the tested strips. Separate the individual parts of the packaging according to material and dispose of them in the correct recycling system.

Observe the local and regional guidelines for disposal.

10. Material composition

Each test strip (5mm x 5mm test patch area) contains the following:

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